

MIRADA

RTX™

Precision confidence for
radiotherapy image workflows




Mirada RTx Imaging Software presents an enterprise solution Image Registration and Visualization package.

RTx is a trusted, accurate and comprehensive image registration and visualization software package for radiotherapy professionals. RTx provides advanced workflows for radiotherapy experts looking for all the benefits of a robust image registration workflow solution without being tied to a costly and complex TPS.

Trusted, comprehensive, image registration and visualization

- Supporting image modalities such as CT, multi-sequence MR, PET and SPECT
- Rigid and modality-specific deformable image registration (DIR)
- Qualitative and quantitative registration quality-check options available



Fully featured for Adaptive

- Adaptive recontouring, dose warping and summation capabilities
- Internal tumour volume auto-contouring based on 4D dataset
- Semi-automatic and manual multi-modality contouring of OAR and target volume

Flexible and customizable solutions

We offer flexible and highly configurable software with various licensing and deployment options. Mirada Medical products are vendor neutral and can integrate seamlessly with your PACs. RTx can be deployed as a standalone workstation or as thin client solution for remote reporting.



Configurable, adaptable, and easy to use

- Send ROIs and other findings to PACS for future recall
- Increase reproducibility and time-saving with Saved Sessions
- Customize the user interface by defining hotkeys, layouts and report templates

Customer service

Mirada Medical is passionate about providing world-class customer service.

Dedicated customer success team delivers support to our customers via telephone support line and email service.

Our experienced clinical and technical specialists will provide you with rapid responses to your queries and work together to swiftly resolve any customer support issue you raise.

Easy Upgrades

Updates are routinely provided at no additional cost to our service customers, meaning you will always have access to our latest products and features.

Training

Available for on-site or remote education, Mirada's experienced team will tailor a training package to suit the needs of your institution.

Available from:
Details for your local representative -

UK Office

Mirada Medical Ltd.
New Barclay House, 234 Botley Road
Oxford, OX2 0HP
United Kingdom
info@mirada-medical.com
+44 (0) 1865 817600

USA Office

Mirada Medical USA, Inc.
1630 Welton Street, Suite 927
Denver, CO 80202
United States
info@mirada-medical.com
+1 303.379.9228